

**RadiantLife** HR<sub>X</sub>  
Hair Rejuvenation Laser Comb



safe

effective

non-drug

non-invasive

scientifically proven  
technology

Technology by  
Greatcore, Inc.,  
241 Polaris Avenue, Mt View, CA 94043 USA  
[www.radiantlasercomb.com](http://www.radiantlasercomb.com)

Low Level Laser Therapy  
for

**THICKER**

**FULLER**



**HEALTHIER  
LOOKING HAIR**


you can use at home



suitable for both young and old, men and women

Distributor

What's all this about?



It's the hair .... of course




**RadiantLife** HR<sub>X</sub>  
Hair Rejuvenation Laser Comb

Low level laser therapy (LLLT) has been used for various ailments for decades in Europe, east Asia, and Russia. More than 2,500 clinical papers have been published with not one reporting any adverse effects. US FDA has approved this technology for joint pain relief and now, hair regrowth.

For the past decade, high-end hair clinics and salons across Europe have been using low level lasers for scalp biostimulation. Hair restoration experts around the world have become confident in LLLT efficacy as an adjunctive tool in the fight against both male and female pattern hair loss. The efficacy has been well established, but it is both costly and inconvenient. Clients are typically charged thousands of dollars for a course of treatment and salon visits are twice weekly.


LLLT for hair application is now brought direct to you in the comfort of your own home at a fraction of the cost with RadiantLife HRx.




RadiantLife HRx uses Low Level Laser Technology for hair rejuvenation by supplying pulsed 655 nm laser light energy to the scalp.

Pulsed laser light energy promotes better cell activities and produces deeper laser light penetration compared to continuous wave.


The laser light is able to penetrate about 8 mm into the scalp tissue, just beyond the hair bulb. The energy is absorbed on a cellular level and photo-biostimulation takes place. Blood circulation to the follicle area is also increased.



Before



After 4 weeks




After 12 weeks

The above subject did not use drugs or any topical lotion during the course of treatment with RadiantLife HRx. Most individuals will experience results within 3 months.

Each treatment session typically takes 15 minutes. 3 to 4 sessions per week is recommended for the first 3 months. Subsequently, 2 sessions per week is advised for maintenance.

RadiantLife HRx can be used to work synergistically with other forms of medically recognized hair treatment.

A 1986 paper on LLLT by the International Laser Group Association states "..... it normalises metabolism of tissues, improves blood cell nutrition, and assures a regular sebaceous secretion .... if we start from a normal quantity of hair bulbs per square cm of the various subjects involved, we can state that there was strengthening of hair bulbs ..... hair growth could be actually ascertained - also it being darker than before".



- ▶ Smaller than palm size, lightweight, and cordless, makes it easy and convenient to use
- ▶ Laser output regulated at optimal power
- ▶ Audible scanning guide and auto off functions
- ▶ Class 3R (3a US ANSI Standard) laser product - suitable for home use

- ▶ The housing is made of tough engineering grade plastics
- ▶ Extremely power efficient with pulsed laser light energy delivery
- ▶ The world's first 1.5V laser comb
- ▶ A single 1.5V AA battery last about 5 months with regular use

